

Kodak Readymatic and Kodak Readypro

Dental Chemicals

When you pair Kodak Readymatic dental chem pack and Kodak Readypro chemicals twin pack with your Kodak dental film, you'll get more than one of the most reliable radiographic combinations for dental practices today. You'll get:

- Sharp radiographic images, works in every processor you use
- Films that clear quickly and chemistry that resists untimely oxidation
- Ready-to-use formulas, eliminating the need for mixing and dilution



The reliability and quality you expect from the most trusted name in imaging



Exceptional solutions for your processing needs

KODAK READYMATIC dental chemicals

Kodak Readymatic dental chemicals are designed for use in all automatic roller-transport dental x-ray film processors. They deliver clear, consistent results for all intraoral and extraoral radiographs. Plus, they:

- Provide excellent film transportability and sensitometry
- Are specifically formulated to limit residue on rollers and in tanks
- Will keep for up to two years, unopened and stored between 40° and 80°F (4.5° to 26.5°C)
- Are available in convenient chem packs

KODAK READYPRO dental chemicals

Kodak Readypro dental developer and fixer are specially formulated to keep rollerless processors running. Plus, they:

- Come packaged in pre-measured bottles; all you have to do is open and pour
- Will keep for up to two years, unopened and stored between 40° and 80°F (4.5° to 26.5°C)
- Are available in handy twin packs

Automated Processing Chemicals	Description	CAT Number
Kodak Readymatic dental chem pack	Two 1-gal bottles of developer and replenisher, two 1-gal bottles of fixer and replenisher	102 8869
Kodak Readymatic dental developer and replenisher	Four 1-gal bottles	877 7534
Kodak Readymatic dental fixer and replenisher	Four 1-gal bottles	871 2879
Kodak Readypro chemicals twin pack	Three twin-packs of 1-qt unit of developer, 1-qt unit of fixer each	117 6262
Kodak RP X-Omat developer starter	To start 8 gal (case of 6 bottles)	133 2642
Kodak RP X-Omat developer and replenisher	Case of 4, 4-liter units	895 6799
Kodak RP X-Omat low odor fixer and replenisher	Case to make 16L (parts A & B)	808 5102
Kodak RP X-Omat low odor fixer and replenisher	Two 5-gal units	886 8804

Manual Processing Chemicals	Description	CAT Number
Kodak GBX developer and fixer twin pack	To make 1 gal developer and 1 gal fixer	190 1859
Kodak GBX developer and replenisher	To make 1 gal	190 0943
Kodak GBX developer and replenisher	To make 5 gal	190 0984
Kodak GBX fixer and replenisher	To make 1 gal	190 1875
Kodak GBX fixer and replenisher	To make 5 gal	190 2485
Kodak rapid access dental developer and fixer	One twin pack of 1-pint unit of developer, 1-pint unit of fixer	183 8374



© Eastman Kodak Company, 2006.
Kodak, Readymatic, Readypro and X-Omat are trademarks of Eastman Kodak Company.

FOR MORE INFORMATION, CALL: 1-800-933-8031 (U.S.A.)
1-800-465-6325, EXT. 35147 (Canada) or visit www.kodak.com/dental
If you are outside of the U.S.A. or Canada, please call your local Kodak company for more information.