


Kodak Ultra-speed

Dental Film

The intraoral film that delivers high-contrast and fine detail.

As a true D-speed film, Kodak Ultra-speed dental film quickly generates the high-quality, detailed images you need to make accurate diagnoses. Ideal for both manual and automated processing methods, Ultra-Speed film is a versatile and easy-to-use solution for your practice.



Enhanced image quality.

Because it accommodates such a wide range of processing and handling conditions, Ultra-speed film is a simple solution for your intraoral imaging needs. Exceptionally tolerant of process variations, this innovative film is designed to render excellent image quality.

The standard in safety.

Ultra-speed film is available with a ClinAsept barrier feature to help reduce the risk of infection. The built-in barrier protects the film packet from contaminants, while its innovative design frees up your time and ensures peace of mind. Simply open the barrier envelope, place the film into the cup and you're ready to start processing.

Kodak Ultra-speed dental film

Available in multiple sizes and formats to meet all of your intraoral imaging needs:

- Available in sizes 0-4.
- Select from flexible, easy-to-open Super Poly-Soft vinyl packets or our sturdier paper products for sizes 0 and 2.
- Bite-wing packets – paper packets with bite-wing tabs attached.
- ClinAsept barrier packets for added safety.
- For optimum results, process your radiographs using Kodak chemicals.

Ordering Information

Size	Description	Code	Packets/Package	REF Number
Size 0	Kodak Super Poly-Soft packets	DF-54	100	122 8840
Size 0	Kodak Super Poly-Soft packets (double film)	DF-53	100	122 8931
Size 0	Kodak Super Poly-Soft packets with ClinAsept barrier	DF-54C	75	144 5360
Size 0	Kodak bite-wing paper packets	DF-48	50	834 8658
Size 1	Kodak paper packets	DF-56	100	127 3747
Size 1	Kodak paper packets (double film)	DF-55	100	127 3721
Size 2	Kodak Super Poly-Soft packets	DF-58	150	165 8194
Size 2	Kodak Super Poly-Soft packets (double film)	DF-57	130	175 3664
Size 2	Kodak paper packets	DF-58	150	149 1737
Size 2	Kodak paper packets (double film)	DF-57	150	149 1752
Size 2	Kodak bite-wing paper packets	DF-40	50	839 3043
Size 2	Kodak Super Poly-Soft packets with ClinAsept barrier	DF-58C	100	171 7131
Size 2	Kodak Super Poly-Soft packets with ClinAsept barrier (double film)	DF-57C	100	898 7877
Size 3	Kodak bite-wing paper packets	DF-42	100	129 6771
Size 4	Kodak occlusal paper packets	DF-50	25	166 6163
Size 4	Kodak occlusal paper packets (double film)	DF-49	25	166 6122

For more information, call **1-800-933-8031** or visit **www.kodakdental.com**



Carestream Health, Inc.

©Carestream Health, Inc., 2009.
Ultra-speed, Super Poly-Soft and ClinAsept are trademarks of Carestream Health.
The Kodak trademark and trade dress are used under license from Kodak.

D5-44 146 7760
2009-07-22